

An Investigation into the Interactions between Various Amyloid-Beta Species and a Formulated Fluorescent Probe

Justin Debrow, Alexandra Hughes, Changqing Chen*

Dept of Chemistry & Physics, Salem State University



Outline

Introduction

Experimental Design

Results

Discussion

Conclusion

References

Acknowledgments



Introduction

- Alzheimer's Disease (AD) is a progressive neurodegenerative disorder that affects approximately 5.7 million US citizens.
- The cause of AD is theorized to be related to the accumulation of various Amyloid-Beta species ($A\beta$) in the brain.
- Curcumin, natural polyphenolic compound, has been shown to have an affinity to $A\beta$, which can serve as a fluorescent probe.
- This research investigated the fluorescence of mixtures of $A\beta$ species and curcumin at different $A\beta$ concentrations.

Asp-Ala-Glu-Phe-Arg-His-Asp-Ser-Gly-Tyr-Glu-Val-His-His-Gln-Lys-Leu-Val-Phe-Phe-Ala-Glu-Asp-Val-Gly-Ser-Asn-Lys-Gly-Ala-Ile-Ile-Gly-Leu-Met-Val-Gly-Gly-Val-Val-Ile-Ala-OH

Figure 1: $A\beta$ 1-42 amino acid sequence with highlighted area of interaction with curcumin





Figure 2: Molecular structure of curcumin, the fluorescent probe



Preparation of $A\beta$ -42 and Curcumin testing solutions:

- 5% curcumin in phosphate-buffered saline (PBS) was used as the fluorescent probe.
- 22.2 μM $A\beta$ -42 aggregates, 22.2 μM $A\beta$ -42 oligomer and 22.2 μM $A\beta$ -42 monomer testing solutions were prepared according to reported procedure.
- Fluorescence measurements were performed using a Fluorescence Spectrophotometer.
- $A\beta$ -42 aggregates images were captured under a Fluorescence Microscope.



Fluorescence Microscopy Results




Figure 3: 10 μL of $A\beta$ Aggregate testing solution with (right) and without (left) 10 μL of curcumin solution captured by a fluorescence microscope

In the absence of the curcumin solution, the $A\beta$ aggregates were indistinguishable from the background.

In the presence of the curcumin, the $A\beta$ aggregates are clearly visible.



Fluorescence of Curcumin in the Presence of 20 μ L A β -42

- At the lowest A β concentration level (20 μ L), all curcumin-A β -42 species mixtures produced relatively similar fluorescence emissions and varied insignificantly from one another.

Fluorescence of Curcumin in the Presence of 20 microL A β -42 Species




Figure 4: Fluorescence of Curcumin in the Presence of 20 μ L A β -42.



Fluorescence of Curcumin in the Presence of 100 μ L A β -42

- Comparing the spectra of the curcumin-monomer and the curcumin-oligomer mixtures, the latter showed a relatively small blue-shift of 15 nm.
- Curcumin-aggregate mixture displayed a much smaller fluorescence intensity.

Fluorescence of Curcumin in the Presence of 100 microL A β -42 Species




Figure 5: Fluorescence of Curcumin in the Presence of 100 μ L A β -42.



Fluorescence of Curcumin in the Presence of 200 μ L A β -42

- Fluorescence of the curcumin-aggregate mixture was quenched. Emission peak was blue-shifted.
- Emission peak of the curcumin-oligomer mixture was blue-shifted and had a slight decrease in peak intensity.




Figure 6: Fluorescence of Curcumin in the Presence of 200 μ L A β -42.



Fluorescence of Curcumin in the Presence of 400 μ L A β -42

- Curcumin in the presence of 400 μ L of A β species showed similar fluorescence spectra to the 200 μ L A β species.





Figure 7: Fluorescence of Curcumin in the Presence of 400 μ L A β -42.




Comparison of Spectra


Fluorescence of Curcumin in the Presence of 20 microL A β -42 Species




Fluorescence of Curcumin in the Presence of 100 microL A β -42 Species



Fluorescence of Curcumin in the Presence of 200 microL A β -42 Species



Fluorescence of Curcumin in the Presence of 400 microL A β -42 Species



Conclusion

- Curcumin interacts differently with soluble monomer, oligomer and insoluble aggregate A β species.
- The fluorescence peak intensities and peak positions changed with the type of A β species as well as their concentrations.
- When there is a greater amount of A β species present, the change in fluorescence intensity and shifting are more evident.
- Under a fluorescence microscope, A β aggregates showed no fluorescence under blue light irradiation. However, upon the addition of curcumin, A β aggregates became visible.



References

- [1] Xueli Zhang, Yanli Tian, Can Zhang, Xiaoyu Tian, Alana W. Ross, Robert D. Moir, Hongbin Sun, Rudolph E. Tanzi, Anna Moore, Chongzhao Ran. "NIRF imaging amyloid and monitoring therapy." Proceedings of the National Academy of Sciences Aug 2015, 112 (31) 9734-9739
- [2] Maresova, P.; Tomsone, S.; Lameski, P.; Madureira, J.; et al. Technological Solutions for Older People with Alzheimer's Disease: Review. Curr. Alzh. Res. 2018, 15.
- [3] Ran, C.; Xu, X.; Raymond, B.R.; Ferrara, B.J.; Neal, K.; Bacskai, B.J.; Medarova, Z.; Moore, A. Design, Synthesis, and Testing of Fluoroboron-Derivatized Curcumins as Near-Infrared Probes for In Vivo Detection of Amyloid- β Deposits. J. Am. Chem. Soc. 2009 (131) 15257-15261
- [4] Alzheimer's Disease: Facts & Figures. <https://www.brightfocus.org/alzheimers/article/alzheimers-disease-facts-figures> (accessed Apr 18, 2020).
- [5] Sou, Keitaro et al. "Loading Of Curcumin Into Macrophages Using Lipid-Based Nanoparticles". International Journal Of Pharmaceutics, vol 352, no. 1-2, 2008, pp. 287-293. Elsevier BV, doi:10.1016/j.ijpharm.2007.10.033.



Acknowledgements

We would like to thank

- Dean Gail Gasparich for providing funding for our research
- Salem State University Chemistry & Physics Department for the continuous support
- Prof. Jason Brown for training on fluorescence microscope
- Prof. Jayashree Ranga for her support

